



# STEP UP! To the Best Stacking Step Stool

The only stool with **ANTIMICROBIAL** protection and designed with **INFECTION CONTROL** prevention as its prime objective:



NO welds and NO hollow legs to trap debris/fluids



Cart Washable and Sterilizable\*



Light Weight Aluminium



1-800-661-5432

[www.StackingStepStool.com](http://www.StackingStepStool.com)



Designed for the O.R.  
Equally Useful Hospital Wide



*Delivering Total Healthcare Value*  
— Since 1935 —


# INTRODUCTION

## Product Description

The IMPERIAL SURGICAL Stackable Stepstools are designed for the O.R., but are equally useful hospital wide. Made of light weight aluminum, our Stackable Stepstools are protected with **an anti-microbial powder coating\***. Furthermore, they can support more weight than most of our competitors (see Figure 1).

- **Improved** model
- **Anti-slip surface** for better stability
- **No rubber corners**: this was a source of previous fluid collection (see the picture on next page, circled in red)
- **No Welds** or large Hollow opening to trap debris or fluids
- **Transport Cart** Available (Figure 2)
- **NEW**: You may now purchase a **fully autoclavable** stepstool without the mat! See last page for order details.

\*Patent Pending



Color-  
code Your  
Stepstools by  
Department!

Can support the weight of half a  
Jeep Wrangler!



Figure 1



Figure 2

**Transport Cart**  
(see last page for  
order details)

# CAN THE COMPETITION COMPARE?

## Antimicrobial Coating


Patent pending technology. These stools are often in contact with moisture, fluid and blood; and will benefit from anti-microbial protection. The benefits include:

- Reduced risks of infections<sup>1</sup>
- Prevents the spread of dangerous micro-organisms and suppresses the growth of mold, bacteria, yeast and mildew
- Provides superior durability and protection for essentially the life of the product.<sup>2</sup>

## No Corner Rubbers

Did you know that the corner rubbers found on some of the competition's stackable stepstools can become a source for bacteria? With time, they erode and create deep cavity channels that can collect fluids and thus, create an environment for bacteria to propagate.

IMPERIAL SURGICAL has removed these rubber corners for your safety.



## Compatible

If you have older models or similar looking competitor stepstools, these newly designed models will stack and inter-lock perfectly.

## High density foot pegs

Our foot pegs are made of denser material than the competition, and will not break or come loose easily. And if they do; we will replace them **free of charge!** Only Transport fees will apply (see Warranty).

## Now Available

Reusable and Single-use (O.R.)\* Anti-Fatigue Mats. They provide comfort, support and reduce leg & back fatigue when standing for extended period.

\*Single-Use Replacement Mats to reduce cross contamination when needed.

## 7 years warranty

We offer one of the best warranty in the market. Our 7 years warranty covers all parts.

## World-Class Customer Service

Your satisfaction is our utmost priority.

1. Anti-Microbial coating provides an added level of protection. However, it should not replace your existing cleaning protocols; it only serves to complement them.

2. Chips, scratches and excessive damage will diminish the anti-microbial properties. If the anti-microbial aspect is of significant importance, we recommend replacement within 7 years or when excessive wear or damage is apparent.

# ORDERING INFORMATION

## Product Details

**Material:** Cast Aluminium

**Weight:** 8.3 lbs (3.8 kg)

**Weight Capacity:** 300 lbs (136 kg)

**Dimensions:** 17.4" x 13" x 4.7" (L x W x H)

**Packaging (1):** 19" x 14" x 6" (10 lbs)

**Packaging (4):** 26.5" x 20" x 15" (45 lbs)



Model Number	Description	Color Mat
OR-36363	Stackable Stepstool with mat	Black 
OR-36363-02	Stackable Stepstool with mat	Dark Blue 
OR-36363-03	Stackable Stepstool with mat	Green 
OR-36363-07	Stackable Stepstool with mat	Orange 
OR-36364	Stackable Stepstool <b>without</b> mat	Autoclavable
Optional Accessory		
OR-3693	Transport Cart	N/A

## NOW AVAILABLE

809-1705	Reusable Anti-fatigue Mat	BX-5
899-1710	Single-Use Anti-fatigue Mat	BX-10
OR-36363-1A	Replacement mat (unit)	
OR-36363-1B	Replacement Foot Pegs	



1-800-661-5432 | 514-631-7988  
 StackingStepStool.com  
 IMPERIAL SURGICAL LTD.